


Best[™] SOFTDISO In-Vivo Dose Verification System

Best SOFTDISO is a stand-alone, user-friendly, workflow driven application that performs In-Vivo dosimetric analysis (IVD for 3DCRT, IMRT and VMAT beams using aSi-EPID portal systems.) It checks the correctness of the treatment by comparing the planned dose (from a treatment planning system) with the applied dose measured from EPID during treatment delivery.

- TEAMBEST GLOBAL

 One Manage of the American Street Street
- Compare dose calculated by TPS and dose reconstructed by the EPID signal
- View all treatment fractions dose comparison tests in the same plot
- View CT scan oriented along beam axis direction
- View EPID profiles aligned below CT scan
- View 2D-Gamma comparison test with gamma % and gamma mean value calculation
- Extract printable report of the tests result
- Display statistics plots and histograms of all the tests performed



Best NOMOS One Best Drive, Pittsburgh, Pennsylvania 15202 USA tel: 412 312 6700 800 70 NOMOS www.nomos.com www.teambest.com